



Dermaforce

ADVANCED RADIOFREQUENCY MICRONEEDLING



PRE AND POST TREATMENT

guidance

3D Dermaforce pre and post treatment care

PRE TREATMENT ADVICE

ONE WEEK PRIOR

- Apply anti-viral cream if prone to cold sores.

3 DAYS BEFORE

- Stop using 'active' skincare (ie. AHA, BHA, Retinol (Vitamin A)).
- Avoid direct sun exposure.
- Excessive facial hair, beards and moustaches need to be clean shaven prior to treatment. Avoid use of aftershaves and perfumes on treatment area.

EXPECTED SIDE EFFECTS*

- Redness · Bruising
- Swelling and dryness

* If you experience any other side effects after your treatment, please contact your clinic therapist.

POST TREATMENT CARE

- Avoid direct sunlight for 24 hours, and no sunbeds during your Dermaforce course.
- Do not apply spray tan or self-tan for 2-3 weeks post treatment.
- Use 3D sunscreen.
- Keep skin dry for 24 hours and avoid scratching
- Avoid active skincare for 3-5 days after treatment.
- No swimming for 24 hours
- No exercise or strenuous activity that will cause sweating.

POST TREATMENT SUGGESTIONS

- Use a clean pillowcase.
- Use hydrating and anti-inflammatory skincare.
- Avoid touching skin directly after treatment.
- Avoid non-mineral makeup for 24 hours.
- Avoid swimming, hot tubs, saunas and active exercise for 48 hours.
- Consider the environment. Hot, humid conditions can contribute to excess inflammation.
- Avoid shaving area for 24-48 hours.

- Avoid heat or steam (ie. cooking/hot shower).
- Do not allow pets to lick the skin.
- Face masks should be disposed of and not reused.
- Do not pick flaky skin.
- Drink plenty of water.



POST TREATMENT INFORMATION

- Up to 24 hours: you will possibly experience erythema and mild swelling.
- 2-3 days: reduced erythema and minor skin dryness and flaking.
- 3-4 days: potential dryness.
- 4-5 days: full recovery.
- 4-6 weeks: improvement in collagen health visible.
- NOTE: 3D Multi-Peptide Serum and 3D Multi-Peptide Cream are recommended for home use.

MAR-DMF-CPPC-REV01



www.3d-aesthetics.co.uk

[f @3DlipoLtd](#) [@3Dlipo](#)



3D-lipo Ltd
Treatment
of The Year



3D-lipo Ltd
Equipment Supplier
of The Year



3D-lipo Ltd
Best UK Based
Manufacturer



MD 704918
ISO 13485:2016 &
EN ISO 13485:2016
Medical Devices
Quality Management